



SPECIMEN CUTTING MACHINE

Specimen Cutting Machine

Model	Image	Description	Price (Baht)
Q-2		Max. cutting Section 35x35mm Speed 2800r/min Grinding Wheel 250x2x32mm Dimension: 395x445x355mm Weight 45kg	26,400
Q-2A		Max. cutting Section 50x50mm Speed 2800r/min Grinding Wheel 300x2x32mm Dimension: 450x515x396mm Weight 67kg	33,700
Q-80Z		Max. cutting diameter 80mm (pipe) Speed 2800r/min Grinding Wheel 350x2x32mm Cutting mode: manual & automatic Display: backlight LCD Dimension: 887x795x600mm Weight 150kg Water tank 60L	153,500
Q-100B		Max. cutting diameter 5-100mm Speed 2100r/min Grinding Wheel 350x2.5x32mm Cutting mode: automatic Display: backlight LCD Dimension: 900x790x600mm Weight 150kg Water tank 60L	173,900
SQ-60		Max. cutting diameter 60mm Speed 2800r/min Grinding Wheel 250x2x32mm Dimension: 690x630x710mm Weight 135kg Water tank 60L	83,700
SQ-80		Max. cutting diameter 80mm Speed 2800r/min Grinding Wheel 250x2x32mm Dimension: 650x715x780mm Weight 119kg Water Tank 60L	91,700



PRICE LIST

SPECIMEN CUTTING MACHINE

Model	Image	Description	Price (Baht)
SQ-100		Max. cutting diameter 100mm Speed 2800r/min Grinding Wheel 250x2x32mm Dimension:680x800x800mm Weight 130kg Water Tank 60L	14,200
LSQ -100		Max. cutting diameter 100mm Speed 2800r/min Grinding Wheel 250x2x32mm Dimension: 680x800x800mm Weight 130kg Water Tank 60L	187,600
LSQ -120		Max. cutting diameter 120mm Speed 2 800r/min Grinding Wheel 350x2.5x32mm Dimension: 690x770x1520mm Weight 330kg Water Tank 60L	204,900
QG - 4		Worktable clamping with quick clamber Max. cutting diameter 65mm Speed 2 800r/min Grinding Wheel 250x2x32mm LED lighting, worktable manual left - right moving, cutting disc manual front -back moving Dimension: 990x 670x1 320mm Weight 2 2 0kg	139,300
QG - 4 A		Rotatable chuck clamping Max. cutting diameter 65mm Speed 2 800r/min Grinding Wheel 250x2x32mm LED lighting, chuck manual rotating, chuck manual left - right moving, cutting disc manual front -back moving Dimension: 990x670x1320mm Weight 240kg	142,800



PRICE LIST

SPECIMEN CUTTING MACHINE

Model	Image	Description	Price (Baht)
LDQ-350		Floor & manual handle type, lighting system Max. cutting diameter 110mm Speed 2300r/min Grinding Wheel 350x3.2x32mm Dimension: 750x1050x1660mm Weight 470kg	346,800
LDQ -350A		Floor & automatic type, Program control, large LCD display, lighting system Max. cutting diameter 110mm Move Capacity of Y -axis 200mm Rotating speed adjustable 0 -2300rpm (0.005 - 10mm/s) Grinding Wheel 350 x3.2x32mm Dimension: 1420x1040x1680mm Weight 500kg	720,800
LDQ -450		Floor & automatic type, Program control, large LCD display, lighting system Max. cutting diameter 160mm Max. moving distance X-Axis 200mm, Y-Axis 300mm, Z -Axis 240mm Rotating speed adjustable 0 -2300rpm (0.005 - 10mm/s) Grinding Wheel 450 x3.5x32mm Dimension: 1420x1040x1680mm Weight 500kg Water Tank 120L	1,690,200
GTQ-5000B		Program control, digital display, Max. cutting diameter 60mm Rotating speed adjustable 500-5000rpm Grinding Wheel 200x0.9x32mm Dimension: 750x860x430mm Weight 126kg	482,500
GTQ-5000A		Program control, digital display, Max. cutting diameter 60mm Rotating speed adjustable 500-5000rpm Grinding Wheel 200x0.9x32mm Dimension: 750x860x430mm Weight 126kg	523,300



BPS INSTRUMENT CO.,LTD.

PRICE LIST

SPECIMEN CUTTING MACHINE

Model	Image	Description	Price (Baht)
DTQ-5		Digital display, Precision cutting, main shaft Rotating speed adjustable 0-600rpm Grinding Wheel 100x0.3x12.7mm Dimension: 415x340x270mm Weight 18.5kg	152900

หมายเหตุ:

1. ราคาดังกล่าวเป็นราคายังไม่รวมภาษีมูลค่าเพิ่ม
2. ราคานี้ฟรีจัดส่งเฉพาะ 0-100 กิโลเมตรแรก กิโลเมตรถัดไปคิดเพิ่มกิโลเมตรละ 30 บาท
สำหรับราคารวมสินค้าตั้งแต่ 100,000 บาทขึ้นไป
3. ในกรณีราคาสินค้าไม่ถึง 100,000 บาท คิด Onsite Service 60 กิโลเมตรแรก 4,500 บาท กิโลเมตรถัดไปเพิ่ม
กิโลเมตรละ 30 บาท
4. สินค้าใหม่ทุกรายการรับประกัน 1 ปี(ยกเว้นอุปกรณ์สิ้นเปลือง)