



Grinding & Polishing Machine

Specimen Grinding & Polishing Machine			
Model	Image	Description	Price (Baht)
P-1		Single disc, 200mm polishing disc, single-speed 1400 rpm, 180W, 220V, 1- Phase, 50Hz Weight : 11kg Packing: 380×490×350mm	15,200
P-2		Floor Type, Double discs, 200mm polishing disc, Single speed 1400 rpm, 180W, 220V, 1-Phase, 50Hz Weight: 47kg Packing : 888×478×930mm	29,500
P-2G		Cabinet Type, Double discs, 200mm polishing disc, Single speed 1400 rpm, 180W, 220V, 1-Phase, 50Hz Weight: 59kg Packing: 888×478×930mm	31,100
P-2T		Double discs, 200mm polishing disc, Sing speed 1400 rpm, 180W, 220V, 1- Phase, 50Hz Weight : 21kg Packing: 700×475×320mm	22,500
M-1		Single disc, 230mm grinding disc, Single speed 700 rpm, 220V, 1-Phase, 50Hz Weight: 13kg Packing: 580×470×520mm	26,400
M-2		Double discs, 230mm grinding disc, Stepless speed Changing 50-1000 rpm, 220V, 1-Phase, 50Hz Weight : 40kg Packing : 380×520×420mm	39,800
MP-1B		Single disc, 200mm polishing disc, Stepless speed Changing 50-1000 rpm, 220V, 1-Phase, 50Hz Weight: 16kg Packing: 380×520×420mm	58,100
MP-1C		Single disc, 200mm polishing disc and 230mm grinding discs, Stepless speed Changing 50-1000 rpm, 220V, 1-Phase, 50Hz Weight: 36kg Packing : 930×530×470mm	71,300



Grinding & Polishing Machine

Model	Image	Description	Price (Baht)
MP-2B		Double disc, 200mm polishing disc and 230mm grinding discs, Stepless speed Changing 50-1000 rpm, 220V, 1-Phase, 50Hz Weight: 40kg Packing: 710×680×330mm	91,800
MP-2		Double disc, 200mm polishing disc and 230mm grinding discs, Double speeds 450&600 rpm, 220V, 1-Phase, 50Hz Weight: 40kg Packing: 710×680×330mm	40,200
MoPao2D		Double motor, Double discs, quick discs shift, 200mm or 250mm, constant speed 300/800 rpm, 380V, 3-Phase, 50Hz Weight : 50kg Packing: 730×725×300mm	81,600
Mopao2DE		Double motor, Double discs, quick discs shift, 200mm or 250mm, Stepless speed changing 60-600, and constant speed 150/300 rpm, 220V, 1-Phase, 50Hz Weight: 55kg Packing: 730×725×300mm	133,200
MoPao260		Double discs, quick discs shift, 200mm or 250mm, constant speed 300/600 rpm, digital display rotating speed, 220V, 1-Phase, 50Hz Weight: 50kg Packing: 710×680×330mm	61,100
MoPao260E		Double discs, quick discs shift, 200mm or 250mm, constant speed 300/600 rpm, digital display rotating speed, 220V, 1-phase, 50Hz Weight: 50kg Packing: 710×680×330mm	101,900
MoPao3S		Single point pneumatic loading, Double discs 300mm, Stepless speed control 5-150 rpm (variable speed). Equipped with automatic grinding & polishing head, 6-hole (30mm) holder, 220V, 1-Phase, 50Hz	367,200



Grinding & Polishing Machine

Model	Image	Description	Price (Baht)
BG-30		Single belt, equipped with cooling system, single speed 1400 rpm, belt size 100×920mm, Surface speed 760m/min, 550W, 380V, 3-Phase, 50Hz Weight: 95kg Packing: 680×630×470mm	75,400
BG-32		Double belt, equipped with cooling system, single speed 1400 rpm, belt size 100×920mm, surface speed 760m/min, 550W, 380V, 3-Phase, 50Hz Weight: 160kg Packing: 690×660×1060mm	104,600
MPTD		Semi-automatic grinding polishing head, rotate speed 50rpm, sample force 0-40N, sample capacity 1-3pcs, 2200V, 1-Phase, 50Hz Weight: 13kg Packing: 183×209×252mm	14,200

หมายเหตุ :

1. ราคาตั้งกล่าวเป็นราคายังไม่รวมภาษีมูลค่าเพิ่ม
2. ราคานี้ฟรีจัดส่งเฉพาะ 0-100 กิโลเมตรแรก กิโลเมตรถัดไปคิดเพิ่มกิโลเมตรละ 30 บาท สำหรับราคารวมสินค้าตั้งแต่ 100,000 บาทขึ้นไป
3. ในกรณีราคาสินค้าไม่ถึง 100,000 บาท คิด Onsite Service 60 กิโลเมตรแรก 4,500 บาท กิโลเมตรถัดไปเพิ่มกิโลเมตรละ 30 บาท
4. สินค้าใหม่ทุกรายการรับประกัน 1 ปี (ยกเว้นอุปกรณ์สิ้นเปลือง)