



Product Code

- UTG-0416 Rotatap Motorized Sieve Shaker,
with 2 Dimensional Motion and Time Adjustment,
for 200 mm (8") Dia. Frame Sieves,
220-240 V 50-60 Hz
- UTG-0416/110 Rotatap Motorized Sieve Shaker,
with 2 Dimensional Motion and Time Adjustment,
for 200 mm (8") Dia. Frame Sieves,
110 V 60 Hz

Standards

ASTM C136

The UTG-0416 Rotatap Motorized Sieve Shaker is designed to be used for 8" diameter sieves. This shaker provides 278 oscillations and 150 taps per minute to produce an effective sieving action. The UTG-0416 provides a 2 dimensional shaking motion and can handle up to 6 X full-height 8" sieves and 13 X half-height 8" sieves. Mounting is not required.



Timer	99 minute, digital
Dimensions	635x710x535 mm
Weight (approx.)	90 kg
Power	1/4 HP