



ZTX-PW7045



Specification		Model
		ZTX-PW7045
Main Body	Sharp stereo erect image over a wide field of view	●
View Head	45° inclined and 360° rotatable binocular head	●
	Both of left and right ocular-tube with diopter adjustment $\pm 5dp$	●
	The interpupillary distance is between 55-75mm	●
Zoom Ratio	1 : 6.4	●
Eyepiece	WF10X/20mm	●
Objective Zoom Range	0.7x-4.5x	●
Magnification	7x-45x	●
Working Distance	100mm	●
Illumination	Natural Illumination	○
	Incident & Transmitted Illumination by 3W LED	●
Power Supply	100V~240V/50Hz	○
Optional Accessory	Auxiliary objective: 0.5x, 1.5x, 2x	○
	Adjustable eyepiece: WF15X/16mm, WF20X/12mm	
	Fluorescent ring lamp: 8W	
	LED ring lamp: 48 or 60 LED bulbs	

